

## Technical Notes

### Technical Note 1 LASER TUNING WITH DIFFRACTION GRATINGS

Light incident on a diffraction grating is dispersed away from the grating surface at an angle which depends on the wavelength, so a grating can be used to select a narrow spectral band from a much wider band. The grating equation,

$$m\lambda = d(\sin \alpha + \sin \beta) \quad (1)$$


(see Figure 1) can be differentiated to give the angular spread (dispersion) of the spectrum:

$$\frac{d\beta}{d\lambda} = \frac{\sin \alpha + \sin \beta}{\lambda \cos \beta} \quad (2)$$

When the grating is operated in the Littrow configuration (in which the light is retro-diffracted; see Figure 2), the equation for the dispersion simplifies to

$$\frac{d\beta}{d\lambda} = \frac{2 \tan \beta}{\lambda} \quad (3)$$


**Figure 1** – The Grating Equation. Here  $\alpha$  is the incidence angle,  $\beta$  is the diffraction angle,  $m$  is the (integral) diffraction order,  $\lambda$  is the wavelength of the light and  $d$  is the spacing between adjacent grooves.




## TECHNIQUES

There are two main methods for selecting a narrow spectral band of light, *Littrow* and *Littman*. These methods are shown in Figures 2 and 3. Littman tuning offers higher angular spread and hence narrower spectral feedback; as the angle of incidence is increased toward 90°, however, the efficiency of the grating drops significantly.

**Figure 2** – Laser tuning using a grating in Littrow mode as the feedback element.



**Figure 3** – Littman tuning using a grating near grazing incidence.



Figures 4 and 5 show typical efficiency curves for two popular plane holographic gratings, with groove frequencies of 1800 g/mm and 2400 g/mm. The curves are measured near Littrow.

Figure 6 gives some insight into the decrease in S-plane efficiency as the incidence angle  $\alpha$  increases.

**Figure 4** – Efficiency curve: 1800 g/mm. Red dashed curve: P-plane; solid black curve: S-plane.




Figure 5 – Efficiency curve: 2400 g/mm. Red dashed curve: P-plane; solid black curve: S-plane.





Figure 6 – Efficiency curve: 1800 g/mm, for several incidence angles  $\alpha$ .



#### FOR FURTHER INFORMATION

For additional information, please [contact us](#).

#### SOME TECHNICAL REFERENCES

- M. Littman, "Single-mode operation of grazing incidence pulsed dye laser," *Optics Letters* **3**, 138 (1978).
- M. de Labacherie and G. Passadat, "Mode-hop suppression of Littrow grating-tuned lasers," *Applied Optics* **32**, 269–274 (1993).
- T. Day, "Continuously Tunable Diode Lasers," *Lasers & Optronics*, 15–17 (June 1993).
- P. Zorabedian *et al.*, "Multiple-prism grating designs tune diode lasers," *Laser Focus World*, 103 (February 1993)
- H. Lotem, "Littrow-mounted diffraction grating cavity," *Applied Optics* **33**, 930–934 (1994).

[back to top](#)

Richardson Gratings  
705 St. Paul Street, Rochester, New York 14605 USA  
Telephone: (585) 248-4100, Fax: (585) 248-4111  
E-mail: [gratings@newport.com](mailto:gratings@newport.com)  
Web Site: <http://www.gratinglab.com>